




iScreen™ Oral

Drug Test Interpretation Quiz



I. MULTIPLE CHOICE

1. What is the length of time required for specimen collection of oral fluids?
- 1 minute
 - 2 minutes
 - 3 minutes
 - 4 minutes

ANSWER: _____

2. The absence of a line next to a drug abbreviation indicates a _____ for that particular drug.
- Non negative
 - Negative
 - Invalid
 - All of the above

ANSWER: _____

3. How long should you wait after a person has eaten or drank before performing the oral specimen collection?
- 10 minutes
 - 15 minutes
 - 1 minute
 - 20 minutes

ANSWER: _____

4. If the result is non negative, how do you secure the specimen for laboratory confirmation?
- Place sponge in a clear plastic bag and send to a lab.
 - Secure cap with security seal and package the device for transport to a lab.
 - Place collection chamber in a clear plastic bag and send to a lab.
 - None of the above.

ANSWER: _____

II. MATCHING (Match the appropriate abbreviation to the drug name)

- | | | |
|---------------------|-------|------|
| 1. Marijuana | _____ | COC |
| 2. Amphetamines | _____ | PCP |
| 3. Cocaine | _____ | mAMP |
| 4. Opiates/Morphine | _____ | AMP |
| 5. Phencyclidine | _____ | THC |
| 6. Methamphetamine | _____ | OPI |

Please fill out your contact information to receive a Certificate of Completion:

Name _____

Company _____

Address _____

Phone (____) _____

III. IDENTIFY THE RESULTS AS EITHER: NON NEGATIVE, NEGATIVE, OR INVALID

